


# Secondary infusions: backpriming technique

Accurate Secondary Mode infusion is dependent upon hanging the Secondary container sufficiently higher than the Primary container.



**Note:** The primary tubing must contain a check valve

1. Prepare secondary container—close the roller clamp on the secondary tubing and spike the secondary container
2. Swab the top of the Y-site port on the primary tubing with appropriate antiseptic and attach the secondary line to the port
3. To backprime, lower the secondary container to a level below the primary container
4. Open the roller clamp on the secondary tubing
5. Allow the fluid to back prime from the primary container into the secondary tubing
6. Close the roller clamp once the secondary tubing is primed and the drip chamber is 2/3 full
7. Hang the primed Secondary container from the IV pole
8. Hang the Primary container LOWER using the hanger(s) provided with the Secondary set
9. Open the vent on the drip chamber if the secondary container is glass or semi-rigid
10. **Pump:** Program the secondary infusion open secondary set roller clamp and start the infusion.  
*When infusion is started, ensure drops are falling in the secondary drip chamber and no drops are falling in the primary drip chamber.*

## Helpful hints

- If flow from the primary container occurs, the primary container may need to be lowered further on an additional hanger.
- Secondary flow rates of 500 mL/h or higher may result in unwanted flow from the primary container.
- Backpriming of glass or semi-rigid containers is not recommended because this may wet the vent.

For product support, please contact Customer Advocacy at **888.812.3266** or email [customerfeedback@carefusion.com](mailto:customerfeedback@carefusion.com).

For technical support, please contact Instrument Technical Support at **888.812.3229**.

For product orders, please contact Customer Order Management at **800.482.4822**.

CareFusion  
San Diego, CA

[carefusion.com](http://carefusion.com)

