



Recommended head height differential

Alaris® System


A. Appropriate head height differential positioning for mini-bags



B. Appropriate head height differential positioning for secondary bottle



C. Appropriate head height differential for intermittent burette delivery



CareFusion recommendations

If necessary, use two hangers to lower primary container to achieve the minimum 9 1/2" of head height differential.

© 2010 CareFusion Corporation or one of its subsidiaries. All rights reserved. Alaris is a trademark or registered trademark of CareFusion Corporation or one of its subsidiaries. IF2084 (1110/500)

CareFusion
San Diego, CA

carefusion.com

