



# DISA DENTAL STUDIO

## CASE PLANNING FOR ALL COSMETIC CASES

(FROM SINGLE ANTERIOR TEETH TO FULLMOUTH RECONSTRUCTION)

To ensure the successful delivery of a case, please make sure that the lab has the following:

### A. DOCTOR'S CHECK LIST: PRE-OPERATIVE

#### 1. DIAGNOSTIC MODEL

- \* P/V for multiple pours; No alginate
- \* No buff or plaster
- \* Die stone and mounting stone only
- \* Impression includes peripheral and free gingival area - (clear) flawless
- \* Opposing poured in die stone, no bubbles, include all posterior teeth
- \* Hammular notches for HIP occlusion, bite sticks, etc. . . To verify canting

#### 2. PATIENT PROFILE

- \* photos - AACD or relax smile and lip line
- \* F Stop, anterior length, 8 & 9 length
- \* Smile design, arch form and contour corrections
- \* Tissue re-contour, gingioplasty - include posteriors



# DISA DENTAL STUDIO

## B. DOCTOR'S CHECK LIST: OPERATIVE

- \* Impressions (check these very carefully)
- \* Opposing model
- \* Pre-op, temp model
- \* Symmetry Bite(Patient standing up)Provide Pictures!!!!!!!
- \* Keep all models
- \* Photos
- \* Instructions, Patient Profile
- \* Photos: pre-op
- \* Temps - feedback. (Work out all problems with the temps and include pictures or models if any changes)
- \* Check temporary for any canting, proportion discrepancies, contour, arch symmetry, length - especially 8 & 9(Provide Pictures)!!!!!!!
- \* Photo of temps 24 hours after prep - lip line and relaxed
- \* Shade, stump shade, mamelon design, surface texture, translucency, etc.
- \* Check clearance with prep matrix
- \* If patient has black triangles or open embrasures, Preps should extend through contact area to Lingual aspect; also check cervical interproximal margin prep.
- \* Elbow preps, not ski jump
- \* Correct prep on Bi's occlusal area, central dissectional groove to buccal
- \* If possible, block out dark dentin areas
- \* Eliminate all sharp line angles
- \* Make clear and accurate records of all comments you feel necessary to assure success with your case